SeQuent[®] Please OTW

- Ⓓ Paclitaxel freisetzender PTA Ballonkatheter
- Ⓖ Paclitaxel releasing PTA balloon catheter
- Ⓕ Cathéter à ballonnet pour PTA à libération de Paclitaxel
- Ⓔ Catéter con balón para ATP con liberación de Paclitaxel
- Ⓘ Catetere a palloncino PTA a rilascio di Paclitaxel
- Ⓒ Paklitaxelavgivande ballongkateter för PTA
- Ⓕ Paklitakselia vapauttava PTA-pallolaajennuskatetri
- Ⓖ Paclitaxel-vrijmakende PTA ballonkatheter
- Ⓕ Cateter balão PTA com libertação de Paclitaxel
- Ⓒ Balónkový PTA katétr s uvolňováním Paclitaxelu
- Ⓒ Balónikový PTA katéter s uvolňovaním Paclitaxelu
- Ⓕ Cewnik balonowy do PTA uwalniający Paclitaxel
- Ⓒ 紫杉醇释放 PTA 球囊导管
- Ⓕ Paklitaxsel salınımı yapan PTA balon kateteri
- Ⓒ Cateter PTA cu balon cu eliberare de Paclitaxel

SeQuent® Please OTW

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SeQuent® Please OTW

SeQuent® Please OTW is available in lengths of 20 mm, 40 mm, 60 mm, 80 mm, 100 mm, 120 mm and 150 mm and in diameters of 1.5 mm, 2.0 mm, 2.5 mm, 3.0 mm, 3.5 mm, 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm and 8.0 mm.

Caution: Before the SeQuent® Please OTW balloon is inflated, the appropriate length and diameter have to be matched with the length and the reference diameter of the target lesion.

Individual treatment

Before using SeQuent® Please OTW, the benefits and risks for each patient must be individually assessed. When establishing the patient exclusion criteria, the risk associated with anti-platelet therapy should be taken into account. Special consideration is required in patients with recent active gastritis or peptic ulcer disease (PUD).

Indications

- De-novo lesions (primary use in the case of stenoses or occlusions) including small vessels
- Restenosis after balloon or stent interventions
- Pre- and post dilatation during peripheral stent implantation
- Acute or impending vascular occlusion

Contraindications

- Intolerance to paclitaxel and/or the delivery matrix (resveratrol)
- Allergy to paclitaxel and/or the delivery matrix (resveratrol)
- Severe allergy to contrast media
- Pregnancy and lactation
- Cardiogenic shock
- Hemorrhagic diathesis or another disorder associated with increased bleeding risk such as gastrointestinal ulceration which restricts the use of platelet aggregation inhibitors and anticoagulation therapy
- Lesions not suitable for interventional treatment
- Vascular reference diameter < 1.5 mm
- Indication for surgical revascularization
- Contraindication for whichever necessary accompanying medication

Possible complications after balloon dilatation

Possible complications may include, but are not limited to:

- Hematoma at the vascular access site
- Pseudoaneurysm
- Arterial perforation or rupture
- Spasm of the treated vessels
- Thrombus formation
- Arterio-venous fistula
- Vascular complications which necessitate a surgical intervention
- Total occlusion of the treated vessel
- Death
- Distal embolism
- Cerebral circulatory disorders
- Restenosis
- Systemic bleeding
- Side effects due to systemic comedication (see corresponding package leaflet)

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- Hypotension
- Infection
- Limb ischemia
- Dissection of the peripheral vessel

Pharmaceutical information on included drug

The very small quantity of paclitaxel in the blood plasma means that known adverse reactions caused by paclitaxel appear less relevant as compared to a systemic administration of paclitaxel. However, yet unknown side effects cannot be ruled out.

Possible undesirable effects

- Allergic or immunological reaction to the drug, similar agents, or to the balloon coating matrix (resveratrol)
- Alopecia
- Anemia
- Gastrointestinal tract dysfunction
- Hematological dyscrasia (including leukocytopenia, neutropenia, thrombocytopenia)
- Abnormal liver enzyme values
- Histological changes in the vascular wall, incl. inflammation, cell damage or necrosis
- Disorders of the cardiac conduction system
- Myalgia/ arthralgia
- Peripheral neuropathy
- Pseudomembranous colitis

Caution:

Paclitaxel has been studied in vitro and in vivo for mutagenic effects with inconclusive results. However, the relevance of the specific potentially genotoxic mechanism regarding a human carcinogenicity risk is currently not known and cannot be completely excluded. This pharmacodynamic action of Paclitaxel acts in interference with the microtubule disassembly. The relevance of this potentially genotoxic specific mechanism for human cancerogenicity risk is currently not known. It is known from the relevant scientific literature that coated, drug-eluting implants (e. g. DES) may cause inflammatory and/or thrombotic reactions. SeQuent® Please OTW does not include any polymer layer, the drug is applied with a degradable matrix on the balloon surface. Although not demonstrated in pre-clinical studies with the SeQuent® Please OTW balloon system, the abovementioned incidences cannot be completely ruled out. Careful analysis of the risk-benefit ratio should be done prior to using the drug-coated balloon.

Allergies

If the patient has a known allergy to any component of the balloon or its coating (paclitaxel or resveratrol) the possibility of allergic/immunological reactions or local inflammatory reactions resulting in possible damage to the vascular wall cannot be excluded.

Interactions with other drugs

The amount of paclitaxel on the balloon surface corresponds approximately to a few hundredth of the quantity usually used in anti-neoplastic treatment, which makes it rather improbable that interactions with other drugs will occur. However, caution should be exercised when concomitantly administering known CYP3A4 and/or CYP2C8 substrates (including terfenadine, cyclosporin, lovastatin, midazolam, ondansetron) or drugs with high plasma protein binding (PPB) (especially sulfonureas, anticoagulants of the coumarin type, salicylic acid, sulfonamides, digitoxin). For

possible interactions of paclitaxel and additional drugs e.g. administered for oncologic indications, the relevant instructions for use should be consulted. A complete study of possible interactions of paclitaxel in combination with comedications has not been established.

Pregnancy

The use of paclitaxel during pregnancy is contraindicated. It is not possible to state the effects of SeQuent® Please OTW on the unborn child. No clinical data are available on the use of SeQuent® Please OTW in pregnant women, reproduction relevant contraindications and risks are unknown.

Instruction for use

Material required to perform a PTA balloon catheter procedure

- PTA balloon catheter for dilation
- Inflation syringe with gauge
- Introducer sheath with dilator
- 0.014 inch, 0.018 inch or 0.035 inch guidewire of suitable length to cross the lesion
- Contrast medium

Prior to use, it is essential that all the instruments, catheters and individual components are carefully checked. An integrity test of the balloon (with the protective balloon sheath left in place), a complete removal of air from the system and a thorough control of the connector seals are mandatory.

Preparation of the balloon catheter

1. The balloon filling fluid is typically a mixture of 50:50 contrast medium and sterile physiologic saline solution or commercially available contrast medium with 200 mg/ml iodine.
2. Draw approximately 3 ml of the contrast medium mixture into a 20 ml or 30 ml syringe. Connect the syringe to the luer adapter of the catheter and hold, with both components pointing downwards. Aspirate. Repeat until no air bubbles can be seen.
3. The syringe system should be filled with a sufficient amount of filling fluid (described in 1.).
4. Remove catheter from dispenser. Remove protective sheath from balloon and transport safety wire from the guide wire lumen. The catheter is now ready for use.

Technique for use

Note: SeQuent® Please OTW is not suitable for coronary use.

Percutaneous introduction techniques and arteriotomy are both suitable when using introduction sets and guide wires.

- Insert an introducer sheath fitted with a hemostatic adapter using standard technique.
- If necessary, the guide wire lumen can be irrigated with physiological saline solution.
- Carefully introduce the guide wire and place the tip distal of the stenosis.
- Predilate the target lesion
- Carefully introduce the fully deflated balloon-tipped catheter over the proximal end of the guide wire and through the hemostatic adapter of the introducer sheath.
- Slide the balloon catheter forward and place the catheter under contrast media injection and x-ray control at the lesion site.
- As soon as the balloon catheter has crossed the stenotic lesion, inflate the balloon to dilate the artery. Depending

on the patient situation and vessel morphology the inflation should be kept (in general) for a period of at least 30 seconds. After dilatation, the balloon must be completely deflated and under x-ray control, pull the balloon tipped catheter back into the introducer sheath. Check the results of dilatation angiographically.

- If necessary, it is possible to use another balloon-tipped catheter via the guide wire that remains in the vessel. In the case of an unsatisfactory primary result, a repeated dilatation of a balloon with larger outer diameter can be done. For prevention of local overdosage it is not indicated to use another SeQuent® Please OTW in the same target lesion.

Use in special populations

The safety and efficacy of the SeQuent® Please OTW balloon catheter to treat supraaortic and cerebrovascular arterial occlusive disease was not examined.

Warnings / Precautionary measures:

- This product should be used only by physicians with experience in angiography and in percutaneous transluminal angioplasty.
- When removing SeQuent® Please OTW from the packaging and when passing the hemostatic valve, great care must be taken to ensure that the balloon system does not become damaged or non-sterile.
- Direct touching of the balloon, wiping the balloon surface or contact with liquids should be strictly avoided, since this may cause delamination of the balloon coating.
- It is possible to rinse the guide wire lumen with sterile/isotonic saline solution.
- Do not inflate the balloon prematurely. The recommended inflation pressure of the balloon must not be exceeded. It is also recommended that a pressure gauge is used to measure the inflation pressure.
- Do not treat the same lesion segment with more than one SeQuent® Please OTW.
- If any resistance becomes apparent at any time during the insertion procedure, the catheter must not be pushed with applied force. The resistance can indicate damage to the balloon catheter. If resistance occurs while advancing through the introducer sheath, the entire delivery system should be pulled back.
- Predilatation with a non-coated balloon is recommended. The catheter should be placed as soon as possible in the target lesion. Extended manipulation of SeQuent® Please OTW can cause delamination of the coating.

Handling and Precautions

The inflated diameter of the balloon should correspond to the lumen of the target artery. Never use a balloon with a larger diameter. Never shake the catheter before use to de-air the balloon. Do not advance the guide wire into the balloon catheter guide wire lumen if resistance is felt without first identifying the cause for the resisting force and taking remedial action. The balloon catheter is 99.9% resistant to the indicated maximum pressure with a confidence interval of 95% (rated burst pressure). Inflation in excess of the maximum inflation pressure is not recommended as the pressure may cause the balloon to rupture, or cause the balloon joints to fail respectively. Should this occur, deflate the balloon and remove it.

Note: Before use, SeQuent® Please OTW must be inspected to ensure that it was not damaged during transport. Do not use if the sterile packaging is damaged or open. Do not re-sterilize.

Note: For treatment to be successful and to ensure long-term satisfactory results, the entire lesion length should be ideally covered by the coated SeQuent® Please OTW balloon section.

Note: Do not bend or squeeze the balloon portion of the SeQuent® Please OTW before the procedure.

Note: Use only the appropriate liquid balloon inflation medium. Failure to use the appropriate fluid (e.g. having a higher concentration of contrast dye) may lead to prolonged in- and deflation times. Never use air or gaseous inflation media.

Note: After dilatation of the vessel segment with the SeQuent® Please OTW stent can be implanted.

Note: The diameter of the balloon should not be less than the reference diameter of the target vessel.

Note: Do not expose the catheter to organic solvents, e. g. alcohol.

Before the angioplasty procedure is started, all equipment to be used for the procedure, including the balloon catheter should be carefully examined to verify its functionality. The appropriate anticoagulant and vasodilator therapy must be provided to the patient prior to the angioplasty procedure and must be maintained during the procedure.

These catheters should only be manipulated in the body while they are under fluoroscopic observation with high quality radiographic imaging equipment. To decrease the potential for the introduction of air into the catheter system, special attention must be paid to ensure that tight catheter connections are maintained. Use the protective balloon sheath until the catheter is being prepared for insertion.

Medication schedule

For all medications please refer to established and current guidelines. On admission to the hospital, the laboratory findings for prothrombin time (PT) and partial thromboplastin time (PTT) as well as the platelet count should be documented for every patient. The following medication schedule is intended as a possible guideline, and should not be seen as a strict medical regimen.

Before balloon inflation:

- Acetylic Salicylic Acid (ASA) 250 to 500 mg 1 x daily p. o.
- Clopidogrel (loading dose) 300 mg/day

During balloon inflation:

- Heparin 70 IU/kg body weight i. v.
- Repeated bolus administration of 2,500 IU heparin i. v., to keep the activated clotting time (ACT) value above 250 seconds

After balloon inflation:

- ASA 75 to 300 mg 1 x daily p. o.
- Clopidogrel 75 mg/d p. o. for 4 weeks (without additional stent) and 12 weeks in case an additional stent is implanted
- Introducer removal normally when ACT value is below 180 seconds and no additional administration of heparin is applied
- In the event of sub-optimal results and/or residual dissections, low-molecular weight heparin can be given subcutaneously 2 x daily for 2 weeks

Note: The stated dosages and instructions must be carefully checked by the treating physician before using the SeQuent® Please OTW balloon catheter. Due to recent developments other procedures which differ from those described above are possible.

Storage

Protect the product from direct sunlight. With proper storage the product may be used until the expiry date indicated on the packaging. The packaged products should not be exposed to temperatures below +10° C and above +40° C.

Disposal instructions

After use, dispose of product and packing in accordance with hospital, administrative and/or local government policy.

Subject to technical modifications.

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