VascuFlex® Multi-LOC

Multiple Stent Delivery System

- Ⓓ Multipler, selbstexpandierender, peripherer Gefäßstent
- Ⓖ Multiple, Self Expandable Peripheral Vascular Stent
- Ⓕ Endoprothèse vasculaire multiple auto-expansible périphérique
- Ⓔ Endoprótesis vascular periférica autoexpandible múltiple
- Ⓘ Stent vascolari periferici multipli autoespandibili
- Ⓔ Ett flertal självexpanderande, perifera kärlstent
- Ⓕ Moniosainen itsestään laajeneva perifeerinen verisuonistentti
- Ⓖ Meervoudig zelfexpanderende, perifere bloedvatstent
- Ⓕ Endoprótese vascular periférica autoexpansível múltipla
- Ⓒ Několicí násobný samoexpandibilní stent do periferních cév
- Ⓒ Mnohohočetny, samoexpandovateľny, periférny cievny stent
- Ⓒ Wieloczołonowy samorozprężalny obwodowy stent naczyniowy
- Ⓖ Мультиплен саморазгъващ се периферен съдов стент
- Ⓖ Többszörös, öntáguló, perifériás ér-stent

VascuFlex® Multi-LOC

DE	Gebrauchsanweisung	3
GB	Instructions for use	8
FR	Mode d'emploi	13
ES	Instrucciones de uso	18
IT	Istruzioni per l'uso	23
SE	Bruksanvisning	28
FI	Käyttöohje	33
NL	Gebruiksaanwijzing	38
PT	Instruções de uso	43
CZ	Návod k použití	48
SK	Návod k použitiu	53
PL	Instrukcja użytkowania	58
BG	Ръководство за употреба	63
HU	Használati utasítás	68

Produktbeschreibung

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VascuFlex® Multi-LOC

Description

The VascuFlex® Multi-LOC system is an innovative stent system for peripheral vascular application that consists of several 13 mm long self-expanding nitinol stents on a single carrier system.

A total of 6 stents are premounted on the coaxial carrier system, and either all or part of these can be placed in the blood vessel as required. To do this, the surgeon releases one or more stents at a specific site (target lesion), moves the entire stent system with the remaining stents up to the next diseased section and there once more releases one or more stents, depending on the situation. This process can be continued as appropriate for the length and number of lesions until all six stents have been implanted.

The VascuFlex® Multi-LOC stent system is introduced into the body through a 6F (2 mm) sheath introducer and via a max. 0.035" (0.89 mm) guide wire and is advanced up to the stenosis that has been treated beforehand by means of a balloon catheter, or up to the occlusion. The self-expanding mechanism of the first stent is triggered by withdrawal of the outer catheter of the carrier system and by the effect of the body temperature. The stent attains its predefined diameter and presses the stenosing material (dissection membranes, calcium deposits or pieces of connective tissue) into the vascular wall by means of the radial force it exerts, thereby restoring the patency of the vascular lumen. If required, additional stents may be implanted in the same vessel, as described above.

The decisive advantage compared to conventional peripheral self-expanding (mono-) stent systems with only one stent is that the surgeon can introduce and release several short stents with only one carrier system. This is advantageous in situations where there are several short high-grade stenoses in an artery, each of which can be covered individually with a short stent, so that the patent segments between the stents remain untouched and the branch-offs of their side branches are not endangered.

The VascuFlex® Multi-LOC system offers a significant advantage when stenting vascular segments subject to strong mechanical strain by bending and stretching, such as for instance the popliteal artery. Here, the stents can be placed one after the other in the case of long stenoses and multiple stenoses (minimum distance of 5 mm recommended), so that even a longer vascular section is covered. Since the individual stents can move freely with respect to each other though, there is nevertheless no risk of stent fractures here, which render the use of long stents in such areas so problematic.

Accordingly, the VascuFlex® Multi-LOC stent system has been developed especially as an alternative to long peripheral

stents, in connection with the use of which such problems occur regularly, entailing the risk of poor morphological and clinical long-term results.

The Stents

Six stents are premounted on a single-hand carrier system with a highly flexible braided outer sheath, thus forming the common VascuFlex® Multi-LOC stent system.

The self-expanding medicut stents of the VascuFlex® Multi-LOC system, which were developed in a so-called "closed cell" design and exhibit only minimal "foreshortening", feature a high degree of flexibility both before release and after implantation in the vessel. The 13 mm long stents are available in four different nominal diameters (5 mm, 6 mm, 7 mm, 8 mm).

Exposed to body temperature, the stent expands to its preshaped nominal diameter after insertion. This high radial force works towards the vessel wall with its deposits and therefore enables the recovery of the patency of the vessel.


Fig. 1 – Diagram of the VascuFlex® Multi-LOC nitinol stent


Fig. 2 – Schematic drawing of the delivery system

Delivery system

The delivery system of the VascuFlex® Multi-LOC system consists of an interior and an exterior catheter. The lumen of the interior catheter is the guide wire lumen. The stents are fitted at the distal end of the delivery system between the interior and the exterior catheter. In order to facilitate the positioning of the stents one radiopaque ringmarker indicates the proximal end of each stent. An atraumatic tip is molded to the distal end of the interior catheter.

The interior catheter is proximally connected to a stainless steel tube. This ends proximally in a Luer-Lock-connector (1) through which the guide wire exits. The exterior catheter is proximally connected to a Luer-Lock slider (2).

The special braided outer sheath of the carrier system ensures, amongst other features, optimal kink-resistance, maximum flexibility, lowest-possible exertion of force for stent implantation and excellent "pushability".

A radiopaque marker band (5) is disposed at the distal end of the outer catheter. The marker band (5) serves to improve the visibility of the outer catheter in the stent area and allows an exact deposition of the stent in the target vessel.

The carrier system is available in two working lengths, namely 80 cm and 130 cm.

Indications

The field of application of the VascuFlex® Multi-LOC nitinol stent is in the treatment of stenoses and occlusions of peripheral arteries.

An important indication for the VascuFlex® Multi-LOC nitinol stent is in the treatment of vascular lesions in segments subject to strong mechanical strain by bending and stretching, such as for instance the popliteal artery. Here, the stents can be placed next to each other in the case of long stenoses and multiple stenoses of such a vessel. Since the individual stents can move freely with respect to each other, the risk of stent fractures is kept to a minimum. The vascular application of the VascuFlex® Multi-LOC nitinol stent is limited to the peripheral vessels including the pelvic arteries. The following indications apply in particular:

Peripheral vessels:

- Insufficient angioplasty result with relevant residual stenosis following PTA
- Reopening of arteriosclerotic vascular occlusions and stenoses, as alternative to treatment with PTA and/or lysis
- Elastic (recoil) stenoses in bypass anastomoses
- Vascular wall dissection

Contraindications

Presently known contraindications for peripheral, vascular

application:

All known contraindications for PTA are also applicable for stent implantation. Contraindications for stent implantation include, but are not limited to:

- Untreated coagulation disorder: coagulopathy
- Existing contraindications for anticoagulant drugs
- Functionally relevant obstruction of the inflow vessel, non-existent outflow or non-existent distal connection
- Stenoses in the vicinity of acute or subacute thrombi
- Highly calcified stenoses that cannot be dilated by means of balloon PTA
- Malposition of guide wire and/or application system
- Verified vascular perforation
- Existing allergies to the materials used, particularly nickel-titanium and tantalum
- Pregnancy (X-rays)

Relative Contraindications

- Renal insufficiency with the risk of renal failure
- Very poor general condition
- Pregnancy

Possible complications of vascular stent implantation

In principle, stent implantations can entail the occurrence of all of the complications known to be associated with surgical procedures, e.g. sepsis, heart failure, kidney failure, coronary ischemia, and death.

- Malimplantation of stents distal or proximal to the lesion
- Inadequate fixation due to selection of too small a stent diameter with possible dislocation
- Acute thrombosis due to lacking perfusion (e.g. spasms or obstructions of other kinds) of the drainage vessels and insufficient anticoagulation
- Vascular rupture with extravascular hemorrhage
- Floating off of atherosclerotic material with peripheral embolism
- Intima injury with dissection
- Formation of an arteriovenous fistula
- Infection-related complications due to failure to maintain sterile conditions
- Surgical intervention to remove the stent rendered necessary
- Hematoma at the puncture site
- Formation of a false aneurysm at the puncture site
- Incorrect positioning technique can result in stent failure when stent is implanted.
- Stent fracture
- Certain complications, e.g. formation of larger pseudoaneurysms, stent dislocation, vascular rupture or

infection, may require a surgical intervention, possibly also under emergency conditions.

- Unwanted overlapping of stents
- Defect of the introduction system

Warnings and Precautions

Please remember that failure to observe the following may result in serious patient complications.

In general

- The stent system should only be used by physicians who are experienced in interventional techniques as PTA or stent procedures or under guidance of physicians experienced in this field. The physicians should be familiar with the complications, the side effects and the dangers that are commonly attributed to stenting.
- The stent system is only suited for single use. Reprocessing and sterilisation of the stent system (or even parts of it) may lead to major faults such as stent system malfunctions. The mechanical properties and material properties can no longer be guaranteed. Risk of infection and risk of contamination may arise.
- The stent system may not be used, if the sterile package has been opened or damaged.
- The stent system may not be used if it is visibly or noticeably damaged.
- The stent system must be used before the expiry date.
- The non-compliance of sterile precautions can lead to infectious complications.
- The VascuFlex® Multi-LOC stent system is not intended for the employment in the heart, in the central circulation system or central nervous system.
- The stents of the VascuFlex® Multi-LOC stent system should not be implanted if other stents from different materials affect the VascuFlex® Multi-LOC stent or are in the vicinity. In this case, electrolytic corrosion and potentially thrombosis formation can be caused.

The delivery system

- The delivery system is not designed for high pressure injections.
- Do not use with Lipiodol or Ethiodol contrast media.
- Do not expose the delivery system to organic solvents (e.g. alcohol).

During the Implantation

- The stent system must always be inserted, moved or retracted over a guide wire (0.035").
- As soon as the release respectively the partial release of the stent has been started it **may not be retracted into the delivery system or repositioned.**
- As soon as the system has been inserted into the vascular system it may only be manipulated under fluoroscopic monitoring. For the X-ray visualization only high-resolution equipment should be used.
- The stent system may not be advanced against strong resistances. Before the procedure can be continued the cause of the resistance must be determined under fluoroscopic monitoring and measures for the retraction

must be considered.

- The blockage of vital lateral branches by the stent is to be avoided.
- A faulty positioning technique can lead to a failure in the release of the stent and thereby to an incorrect implantation.

The Stent

- After removal of the safety pin (4) the stents should be easy to release; i.e. the hand wheel (3) should be easy to turn. Under no circumstance try to release the stents by exercising increased force and make sure that the safety pin (4) has been removed. In case of malfunction remove the complete system and replace it by a new system.
- As soon as a stent is partially released it may not be retracted into the delivery system. Also it may not be repositioned.
- If several lesions are to be treated during one procedure, first the distal lesion and then the proximal lesion should be stented. Stenting in this order avoids having to pass already positioned stents again, thereby possibly dislocating them.
- The renewed passage of a previously implanted stent should only be exercised if absolutely necessary and then with great care.
- In order to prevent an overirritation of the vessel wall and to maintain the flexibility of the vessel it is best advice to position the stents not end to end and not overlapping.
- If the size of the stent is not matched exactly to the native vessel the stent might migrate if the diameter has been chosen too small due to the insufficient radial force. In the case of a diameter chosen too large the intima can be damaged due to too high radial force. The stent placement is to be chosen so as to completely cover the hemodynamically relevant area of the lesion (dissection, flap, or recoils). For the correct choice of stent diameter please confer the selection table.

Implantation Procedure

Artificial stents must always be implanted by surgeons that are trained and familiar with the specific techniques of such operations. Please also observe the following points when implanting the VascuFlex® Multi-LOC stents:

Preparation

- The percutaneous placement of the stent in a stenotic or obstructed vessel is to be exercised under fluoroscopic monitoring in an angiographic cath-lab.
- The preparation of the patient as well as the sterile measures to be taken are to be the same as for any angioplastic catheter treatment.
- The accompanying medication is determined by the treating physician, including the necessary anticoagulation required during and after the treatment. These measures must be aligned to the most up-to-date scientific knowledge and must be adjusted by the treating physician to the respective patient and the specific situation.

Angiography

- The extent of the lesion(s) and the collateral circulation

must be established employing angiography .

- If there is a thrombus formation or a suspected thrombus formation a thrombolysis should be exercised using a clinically customary and proven standard method.
- The access vessels must be sufficiently patent or re-channelled, before the intervention can be continued.
- The lesion should be identified and marked under fluoroscopy and under consideration of the distal part of the stenosis or stricture.

Pre-dilation of the Lesion

- It is advised to pre-dilate the lesions intended to treat by standard PTA technique in order to be able to determine the correct size of the stent.
- The balloon catheter is removed and the access to the lesions is maintained by the remaining guide wire.

Selection and preparation of the stent system and compatibility with the accessories

1. Selection of stent size

- **Please note:** Should more than one stent be necessary to cover the lesion the distal stent should be implanted first, followed by the stent(s) to be placed proximally.
- In order to avoid any irritation of the vessel wall and to maintain the flexibility of the vessel it is best advice to position the stents not end to end and not overlapping. A minimum distance of 5 mm between two stents is recommended to change as little as possible the biomechanical properties of the artery.
- If the dilation of a vessel is necessary prior to placement of the stent from the medical point of view, then the balloon diameter must not be larger than the nominal diameter of the stent.
- In order to determine the correct diameter of the stent, the diameter of the reference vessel (proximally and distally to the lesion) must be measured. It is necessary to select a stent which expanded nominal diameter is not more than one millimeter bigger than the largest diameter of the reference vessel so that the stent does not overdilate the vessel, which could lead to intima proliferation with stenosis formation.
- For the correct selection please take note of the following table 1:

Table 1 – Stent-Selection Table

Reference Vessel diameter	Nominal diameter of the stent in expanded status
4.0 - 5.0 mm	5 mm
5.0 - 6.0 mm	6 mm
6.0 - 7.0 mm	7 mm
7.0 - 8.0 mm	8 mm

Usable length of the delivery system: 80 or 130 cm
 Recommended introducer sheath: min. 6F (2 mm)
 Recommended guide wire: max. 0.035" (0,89 mm)

2. Preparation of the delivery system

- First of all please make sure that the sterile package has not been damaged or opened.
- Open the sterile pouch and remove the blister under sterile conditions.
- Take out the delivery system from the blister. Check it for visual or noticeable damages. Check whether the stents are mounted properly on the delivery system.
- **Please note:** If there is reasonable doubt that the sterility or the functionality is impaired, the system may not be used.
- Before the use of the delivery system it must be completely rinsed with heparinized saline solution through both proximal accesses.
- When rinsing via the Luer-Lock slide valve (2), so much heparinized saline solution must be injected that the rinsing liquid emerges visibly at the tip. Only if this is the case the lubrication between the sheath introducer and the filling material and easy placement of the stents is ensured.
- Then the central guide wire lumen is rinsed by connecting the syringe to the proximal Luer-Lock connector (1).

Stentimplantation

1. Insertion of the delivery system

- Using common techniques, such as the Seldinger technique, a percutaneous access is created. Then an introducer sheath compatible with the delivery system (min. 6F / 2 mm) is inserted.
- Through the introducer sheath a max. 0.035" (0,89 mm) guide wire of suitable length is advanced up to and beyond the stenosis/stricture.
- First of all a balloon catheter is inserted over the guide wire and the lesion to be treated is pre-dilated. Then the balloon catheter is removed. The guide wire remains in place to maintain access to the lesion.
- **Please note:** All manipulations with the delivery system or other devices must always be exercised over the guide wire.
- Always use an introducer sheath with a haemostatic valve. For contra-lateral accesses it is advised to use a long introducer sheath specially designed for contra-lateral accesses which simplifies the crossing of the bifurcation.
- In order to remove air bubbles from the guide wire lumen and the interior catheter these are once again rinsed with heparinized saline solution.
- The delivery system is advanced to the lesion over the guide wire under permanent fluoroscopic control. The exact positioning of the stent in the stenosis/stricture is performed with the aid of the radiopaque markings on the proximal side of the stent as well as the marker band (5) at the tip of the outer catheter.
- **Warning:** The stent system must not be advanced against strong resistance. Before the procedure can be continued the reason for the resistance must be determined under fluoroscopic monitoring, and measures for the retraction must be considered.

2. Release of the stent

- Before release of the first stent, the red safety pin (4) at the handle of the release system must be pulled out. Subsequently the exact position of the stent in the stenosis / stricture should be verified. The hand wheel (3) shall not be turned before the safety pin (4) has been removed.
- **Warning:** The stents should be easy to release by turning the hand wheel (3). Interrupt the procedure if an unusual effort is necessary. A change of system is necessary.

After stent implantation

- After the complete release of the intended number of stents the complete delivery system is removed slowly and carefully under fluoroscopic control over the guide wire out of the body. Thereby the guide wire remains still in position.
- **Please note:** Stents that are still in the carrier system may only be placed in the course of the ongoing procedure and in the same patient. If not all stents are required, then the carrier system along with the rest of the stents still remaining in the sheath introducer must be disposed of. Reprocessing / re-sterilization of the carrier system with the non-required stents is not permitted.
- **Please note:** Should the tip of the delivery system get caught in a distal stent segment when pulled back it is to be freed from it by careful movement.
- By conventional execution of a post angiography it is checked whether all implanted stents have completely expanded.
- Has a stent not completely expanded at any position within the lesion, an additional balloon dilation (PTA technique) can be exercised. For this a balloon catheter of suitable size is chosen and the lesion dilated in common technique. The inflated diameter of the balloon has to correspond with the reference vessel.
- **Caution:** Never expand a stent with a balloon which is larger in diameter than the expanded diameter of the stent.
- Has the procedure been successful, the guide wire and the introducer sheath can also be removed and the haemostasis of the access site can be induced.
- Should the result of the stenting not yet be satisfactory, further measures to be taken must be decided.

MR – Compatibility

The stents are MR-compatible. They will not be influenced by MR and cause only low or none image artifacts. Nevertheless, it is recommended to carry out a MR examination only about six weeks after the implantation to avoid a dislocation. The current standard of knowledge regarding magnetizing force and adjustment to the magnetic field should be considered to receive optimal results.

Storage and packaging

VascuFlex® Multi-LOC is a single-use product that has been sterilized with ethylene oxide. It is supplied in sterile packaging and will remain sterile as long as the package is not opened or damaged. The product is free of pyrogen.

The stent system must be kept in a dark, dry and cool place. Direct sun exposure must be avoided – heat could have a considerable, negative influence on the functionality of the product.

Do not expose the stent system to ionizing radiation or ultra-violet light.

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